

HEPARIN TUBE (green top)

USES:

LFTS
RFTS
ELECTROLYTES
SERUM LDH
TOXICOLOGY
SERUM PHOSPHORUS
SERUM MAGNESIUM
CRP
DIRECT/INDIRECT BILIRUBIN
HBSAG/ANTI HIV/ANTI HCV
HEP A E ANTIBODIES
DENGUE NS1
URIC ACID
AMYLASE/LIPASE
RA FACTOR
CPK
INDIRECT COMB TEST
BLOOD GLUCOSE
OSMOTIC FRAGILITY TEST



Fb/Nurse-Info

IDZZZZZZZZZZZZ

EDTA TUBE (PURPLE TOP)

USES:

CBC

BLOOD SMEAR

ESR

RETIC COUNT

MALARIAL PARACITE

HBA1C

DIRECT COMBS TEST

G6PD

HB ELECTROPHORESIS

RED CELL FOLATE LEVELS

MONOSPOT TEST

BLOOD GROUPING

CROSS MATCH



Fb/Nurse-Info

JDzzzzzzzzzzzz

CITRATETUBE (BLUE TOP)

USES:

PT/INR

APTT

D DIMERS

FDPs



JDzzzzzzzzzzzz

Fb/Nurse-Info

GEL TUBE (YELLOW TOP)

USES:

CARDIAC MARKERS

FERRITIN

TUMOR MARKERS

SERUM CORTISOL

ANTIBODIES ANA

ALPHA FETO PROTEIN

CEA, CA125, PSA, CA19

HCV PCR

VIT D

SERUM OSMOLALITY

TFTS



Fb/Nurse-Info